


Timion Child Equipment

LONG SEATER

3. Adjust the angle if a child cannot handle the pull on the hamstrings: Slide the backrest upwards and put a bolt under the seat at the desired height.

2. Adjust knee straps

3. Adjust foot plates



Screws for angle adjustment are stored here.

*If you are unable to assemble the long seater for any reason or if you are unable to deliver it to the child's house, please **do not hesitate to call Timion. We are happy to assist.***

4. Position the child on the seat. Fasten the lap strap



PO Box 4404, Aston Bay 6332, South Africa
Tel: 078 456 1512/ 1870 info@timion.org
Fax: 086 561 7173 www.timion.org

This child is precious, of considerable value, held in honour, and especially dear! (= Gr. timion)

Elbe - 0837504431
Daniel - 0784561512
danika.meyer@gmail.com